



CIB Genomics Core Lab
Blue Pippin size selection

Date Received: _____

Date _____
 Project Name _____
 Contact _____
 Department _____
 UMC _____
 OCC # _____
 E-mail _____
 Phone _____ Fax _____

Ninglin Yin, Genomics Group Leader
Center for Integrated BioSystems
 4700 Old Main Hill
 Logan, UT 84322-4700
 (435)797-8218 or (435)760-3688
 Fax: (435)797-2766
 e-mail: usudna@aggiemail.usu.edu
 Web site: <http://www.biosystems.usu.edu>
Billing Questions: Sue McCormick (435) 797-7510

Please sign/date acknowledgement below:
While the CIB strives to maintain high quality results, if results are not found to be acceptable, investigator will still be billed for all runs and reruns of samples. Lab fees may be negotiated.

Principal Investigator:
 Name (print): _____
 (Sign): _____ Date: _____

Important points regarding samples

- The maximum volume that can be loaded is 30 µl. Samples delivered in volumes larger than 30 µl will be divided between multiple lanes and charged accordingly.
- Samples should be in a low salt buffer with any DNA-binding proteins/enzymes removed. Please clean up enzyme digests, PCR products, etc. first.

Please check box to have the samples analyzed on the TapeStation after size selection. Additional charge is \$5 per sample plus \$5 per TapeStation run.

Charge is \$15/sample.

							Lab Use Only		
	Sample Name	Type of DNA (PCR product, seq library, etc.)	buffer	Conc. (ng/µl)	Total Volume	Size range desired			
1									
2									
3									
4									
5									
6									
7									
8									

_____ **Estimated Price for Project**